

Progress in the development of capillary Atmospheric Pressure Ionization Methods (cAPI)



Physical & Theoretical Chemistry

Wuppertal, Germany

Institute for Pure and Applied Mass Spectrometry

Sonja Klee; Valerie Derpmann; Hendrik Kersten; Walter Wißdorf; David Müller; Yessica Brachthäuser; Sebastian Klopotoski; Klaus Brockmann; Thorsten Benter

Introduction

Challenges

- Design AP ion sources which generate ion populations closely representing the neutral analyte composition *in the source*
→ Investigations of *ionization mechanisms*
- Design AP ion sources which operate kinetically controlled
→ Ionization within *reactive matrices* (e.g. effluent from smog chambers, ambient air)
- Design AP ion transfer stages which allow carefully controlled transport of ions from the source to the analyzer region
→ *Optimization* of ion transfer efficiency (instrument sensitivity)
→ *Controlled* ion activation (reagent ion reactivity)

State of Knowledge:


- All *common* API sources are generating ion populations, which are thermodynamically controlled, *i.e.*, all ion concentrations present are *equilibrium* concentrations [1].
- Thermodynamic control is governing the thermal ion source chemistry since the ion residence time in the source enclosure is generally ≥ 10 ms. Most bi- and termolecular ion molecule reaction rates are collision controlled leading to half lives of $< 10^{-5}$ s (first order rates $> 10^5$ s $^{-1}$) when neutral reactants (*e.g.*, H₂O, dopants, solvents) are present at ppmV mixing ratios.
- Thermal equilibrium ion concentrations do rarely reflect the true neutral analyte/matrix composition. For example most low mass ions are generally present in the source as charge bound clusters [2].
- In the context of the present framework, “reactive” neutral matrices are depicted as containing ppmV (or above) mixing ratios of either *polar* compounds such as water or LC solvents (CH₃OH, CH₃CN) as well as molecular oxygen, due to its bi-radical character (\dot{O}_2)

Methods

Experimental Setup

- MS:** Bruker esquire6000 ion traps, Bruker micrOTOFs, TOFWerk HR-ToF with custom inlet and pumping stage, Thermo Exactive Orbitrap
- Ion Sources:** Custom APPI, APCI, APECI, and APLI stages. Laminar flow ion source (LFIS). Bruker Apollo and MPIS [3] sources for comparison studies
- Capillaries:** Borosilicate glass, stainless steel, aluminum
- Gases:** Nitrogen (5.0), synth. air (5.0) < 1 ppmV H₂O

General Considerations




Conclusions

- Capillary ionization methods are currently operated in the gas phase, *i.e.*, require *vaporization* of the neutral analyte before ionization occurs
- For *direct cAPPI* (no dopants present) kinetic control is achieved. Thus complex mixtures with highly reactive neutrals present can be accurately analyzed. This is of particular interest to applications in atmospheric/environmental chemistry.
- The same arguments apply for cAPECI due to the high reactivity of superoxide
- cAPCI in positive mode using liquid point electrode configurations is fully compatible with cAPI – this method benefits strongly from a dedicated ion activation stage (within domains 4-6), since the reagent ions are fully thermally equilibrated [H+(H₂O)_n]⁺ clusters.

Thermodynamic vs. Kinetic Control

If the residence time of ions in domains 1-5 (cf. Figure 1) is significantly exceeding the half lives of one or several ion/molecule reactions, then $\Delta_R G$ governs the ion concentration distribution.



Kinetic control
The reaction system is *not* allowed to equilibrate before mass analysis.

The ion population resembles closely the neutral analyte composition in the matrix.

$A_g^+ + B_g \rightleftharpoons C_g^+ + D_g$

Thermodynamic control
The reaction system has time to equilibrate, driven by $\Delta_R G$. The ion population resembles the energetically lowest ion system. This may represent the neutral analyte composition, see insert in Fig. 2, but this is rarely the case.


$A_g^+ + B_g \rightleftharpoons C_g^+ + D_g$

Typical Examples

- Radical cations M^{•+} are not oxidized in the presence of O₃ or OH to yield *e.g.* [M+O₃]^{•+}
- Aromatic hydrocarbons, *e.g.*, toluene, do *not* act as proton donors in complex matrices (and are thus *not* “suppressed”)
- High IP analytes do *not* transfer charge

Typical Examples

- H₂O⁺ ions are rapidly converted to yield a thermally equilibrated distribution of [H+(H₂O)_n]⁺ clusters
- Low mass analyte ions are generally clustering, mostly with background matrix water or polar LC solvents
- Ion “suppression” phenomena
- “Matrix effects”



cAPPI (see MOH PI 3:30 pm and MP15#276)


Purpose 1: Minimize the contact time between ions generated at high pressure and reactive neutrals such as radicals present in the gas matrix ($t < 200$ μ s) thus pushing the chemical system towards kinetic control.

Approach: Miniature VUV lamp mounted windowless onto the inlet capillary [4].

Design:

Results:

- Mixture of p-xylene (D) and a 1,4-unsaturated dicarbonyl (A): Progress of the protonation reaction $D^+ + A \rightarrow D-H + [A+H]^+$ for typical APPI dwell times (> 5 ms; bottom left), in contrast to the kinetically controlled ion distribution obtained with cAPPI (< 0.13 ms; top left) [4].
- Pyrene radical cations (m/z 202) are oxidized in the presence of OH-radicals using APPI, but hardly with cAPPI (right) [5]



Purpose 2: Highly efficient direct (*no dopant*) VUV light driven generation of analyte radical cations in an ultra-clean environment. Application to GC-cAPPI-Orbitrap MS.

Approach: GC capillary effluent → transfer line → heated capillary ducts → APPI ionization chamber → capillary ducts → Orbitrap inlet capillary

Design:

Results: Low fg on-column sensitivity for various analytes (as M^{•+}) See MOH Photoionization talk 3:30 pm / Room 103

Capillary Ionization Methods


cAPCI (see MP15#274/5 and #283)

Purpose: Establish a selective *protonation* method for GC hyphenation without the adverse effects accompanied by the interaction of the hot Corona plasma zone with the analyte gas stream [2].

Approach: Thermally equilibrated proton bound clusters [H+(H₂O)_n]⁺ are generated as reagent ions upstream of the capillary duct. The corona is operated with a liquid point electrode (highly purified, optionally acidified water pumped through a quartz capillary at nano-flow rates) and a funnel shaped plane counter electrode merging at the exit with the inlet capillary duct. Since the reagent ion clusters ($n = 5 \dots 9$) exhibit rather low proton acidities, an ion activation stage (domains 4-6) is required for high sensitivities towards less basic compounds (see poster MP15#283).

Design:

Results: cAPCI is a protonation-only API method. Charge transfer as well as oxidation reactions are generally not observed. The figure below shows the TIC for a GC-cAPCI-TOF MS analysis of a fatty acid ester mix (left panel). The corresponding mass spectra consist solely of a strongly concentration-dependent distribution of protonated molecule, proton bound analyte water clusters, and proton bound dimer ions, as expected (right panel). Carefully controlled ion activation in a dedicated rf funnel device (see poster M15#283) leads to exclusive [M+H]⁺ formation.




CAPECI (see MP15#277/8)

Purpose: Establish a selective kinetically controlled ionization method for analytes with elevated electron affinity or proton acidity [6]. Samples are either delivered from a GC with transferline stage or are directly sampled from ambient air. Future applications will include spatially resolved desorption stages coupled directly to the capillary duct (see below).

Approach: Selective generation of thermal electrons via the photo electric effect with long-term stability using anodize aluminum as photo emissive material. Quantitative capture of the photo electrons by oxygen to yield superoxide, O₂⁻, as reagent ions.

Design:

Results: Typical minimum analyte/reagent reaction times in conventional APPI sources are > 10 ms. As shown below for 2,4-dinitro-toluene even a reaction time of 5 ms still leads to the formation of the thermodynamically preferred deprotonated and further oxidized molecule (left panel). With reaction times of < 1 ms, only the molecular radical anion is detected (right panel).



References

- Dzidic, I.; Carroll, D.J.; Stillwell, R. N.; Horning, E.C. “Comparison of Positive Ions Formed in Nickel-63 and Corona Discharge Ion Sources Using Nitrogen, Argon, Isobutane, Ammonia and Nitric Oxide as Reagents in Atmospheric Pressure Ionization Mass Spectrometry” *Anal. Chem.* 1976, 48, 1763-1768.
- Klee, S.; Albrecht, S.; Derpmann, V.; Kersten, H.; Benter, Th. „Generation of Ion Bound Solvent Clusters as Reactant Ions in Dopant Assisted APPI and APLI” accepted for publication as “paper in forefront” in *Anal. Bioanal. Chem.* 2013.
- Schieweck, R.; Lorenz, M.; Giese, R.; Brockmann, K.J.; Benter, Th.; Gäb, S.; Schmitz, O.J. “Development of a multi-purpose ion source for LC-MS and GC-API MS” *Anal. Bioanal. Chem.* 2008, 392, 87-96.
- Kersten, H. “Development of an Atmospheric Pressure Ionization source for in situ monitoring of degradation products of atmospherically relevant volatile organic compounds” Dissertation, Faculty of Mathematics and Natural Sciences, University of Wuppertal, Germany, 2011. <http://elpub.bib.uni-wuppertal.de/eodocs/dokumente/fbc/chemie/diss2011/kersten>
- Kersten, H.; Funcke, V.; Lorenz, M.; Brockmann, K.J.; Benter, Th.; O'Brien, R. “Evidence of Neutral Radical Induced Analyte Ion Transformations in APPI and Near-VUV APLI” *J. Am. Soc. Mass Spectrom.* 2009, 20, 1868-1880.
- Derpmann, V.; Albrecht, S.; Benter, Th. “The role of ion-bound cluster formation in negative ion mass spectrometry” *Rapid Commun. Mass Spectrom.* 2012, 26, 1923-1933.

Acknowledgement

Financial support is gratefully acknowledged

- Bruker Daltonics (Bremen, Germany)
- Thermo Fisher (Bremen, Germany)
- iGenTraX UG (Haan, Germany)
- DFG – German Research Foundation
- DAAD – German Academic Exchange Service
- BUW – University of Wuppertal, Germany
- iPAMS – Institute for Pure and Applied Mass Spectrometry, University of Wuppertal, Germany